



Digital Titrator-patron med kvicksilvernitratt, 0.2256 N

Produktnr:

1439301

SEK Pris:

Kontakta oss

Tillgänglig

Farligt



Bruksfärdig titreringslösning för bestämning av klorid (kvicksilvernitrattmetoden), 10-8000 mg/l Cl⁻.

For chloride determination by Digital Titrator titration with mercuric nitrate, using diphenylcarbazone indicator. Hach method 8206. Digits = mg/L Cl⁻ in 100 mL sample. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Mätområde:	10 - 160 mg/L Cl ⁻
Metodnamn:	Mercuric Nitrate
Metodnummer:	8206
Parameter:	Chloride
Plattform:	Digital Titrator
Titrant:	0.2256 N Hg(NO ₃) ₂