



Kuddar med pulver för syrareagens, 25 mL, pk/1000

Produktnr:

104228

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



For determination of high range silica by the silicomolybdate method (Hach method 8185) Range: 1-100.0 mg/L SiO₂. Also used in the iodometric method (Hach method 10100) for the determination of high range hypchlorite (50-150 g/L as Cl₂). Pack of 1000 powder pillows.

Specifikationer

Metodnamn:

Silicomolybdate

Plattform:

Powder Pillows