



MolyVer® 6 molybdenum reagent AccuVac® ampules, 25/pk

Produktnr: 2522025

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Föreslagna ersättningar:

- 2604100



For determination of HR Molybdate Molybdenum by the Mercaptoacetic Acid Method. Hach method 8036. Range: 0.2 - 40.0 mg/L Mo. Pack of 25 ampules. Also requires CDTA Solution, 0.4 M, Cat. No. 26154-36.

Fast

AccuVac® Ampules provide quick results. Just snap, fill and read - saving you valuable time.

Easy

No need to measure sample or reagents. The evacuated ampules draw in the correct sample volume, so anyone can obtain reliable results with minimal effort.

Accurate

Optical quality evacuated vials assure precise measurements and reagent stability.

Use them anywhere you need results.

Ampules are both reagent containers and photometric measurement cells. No need to transport reagents or samples.

Specifikationer

Antal tester :	25
Mätområde:	0.2 - 40.0 mg/L Mo
Metodnamn:	Mercaptoacetic Acid
Parameter:	Molybdenum
Plattform:	Accuvac