



Reagent, hydrazine, AccuVac® ampules, pk/25

Produktnr: 2524025

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



Föreslagna ersättningar:

- 179032

For determination of Hydrazine by the p-Dimethylaminobenzaldehyde method, using AccuVac® Ampules. Hach Method 8141. Range: 4 - 600 µg/L. Pack of 25 ampules.

Fast

AccuVac® Ampules provide quick results. Just snap, fill and read - saving you valuable time.

Easy

No need to measure sample or reagents. The evacuated ampules draw in the correct sample volume, so anyone can obtain reliable results with minimal effort.

Accurate

Optical quality evacuated vials assure precise measurements and reagent stability.

Use them anywhere you need results.

Ampules are both reagent containers and photometric measurement cells. No need to transport reagents or samples.

Specifikationer

Antal tester :	25
Instrument:	All except DR/850, DR/820, PC II
Mätområde:	4 - 600 µg/L N ₂ H ₄
Metodnamn:	p-Dimethylaminobenzaldehyde
Metodnummer:	8141
Parameter:	Hydrazine
Plattform:	Accuvac