



Karusellhållare 7 x 1 cm

Produktnr: LZV902.99.00002

SEK Pris: Kontakta oss

Skickas inom 1 vecka

The carousel sample changer allows up to seven sequential measurements for mini series or enzymatic methods.

Management of blank sample and analysis sample for multiple cells.

The carousel program wizzard allows different sequences for blank and analysis samples.

Automatic result calculation for time depending absorbance readings in a unique sample.

The free programming mode allows the calculation of concentration depending on time-triggered absorbance readings.