



Digital titratorpatron, natriumhydroxid, 8,00 N

Produktnr:

1438101

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For determination of acid concentration by Digital Titrator titration. Digits x 0.01 = meq

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Mätområde:	10 - 4000 meq/L Acid
Metodnamn:	Phenolphthalein End Point
Metodnummer:	8233 (Base)
Parameter:	Acid-Base
Plattform:	Digital Titrator
Titrant:	8 N NaOH