



CDTA Digital titratorpatron, 0,800 M

Produktnr:

1440301

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Cyclohexanediamine tetraacetic acid standard solution. For Calcium and Total Hardness determination by Digital Titrator titration. 0.800 M. Digits = mg/L CaCO₃ in 100 mL. Digits x 0.001 = mmole.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Metodnamn: CDTA Titration

Plattform: Digital Titrator