



Digital titrator cartridge, sodium hydroxide, 0.9274 N

Produktnr:

1484201

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



For volatile acids determination by Digital Titration, 100 - 2400 mg/L Acetic Acid

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Mätområde:	100 - 2400 mg/L CH ₃ COOH
Metodnamn:	Phenolphthalein End Point
Parameter:	Volatile Acids
Plattform:	Digital Titrator
Titrant:	0.9274 N NaOH