



Digital titrator cartridge, phenylarsine oxide (PAO), 0.00564 N

Produktnr: 199901

Farligt

UTGÅNGEN ARTIKEL
Den här artikeln är inte längre tillgänglig.



For low level chlorine determination by amperometric titration, using the Digital Titrator

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

EPA compliant:	Ja
Mätområde:	0 - 1000 µg/L Cl ₂
Metodnamn:	Amperometric, Forward
Metodnummer:	10026
Parameter:	Chlorine, Free
Plattform:	Digital Titrator
Titrant:	0.00564 N PAO