



Free and total chlorine test strips, 0 - 10 mg/L, 50 tests

Produktnr: 2745050
SEK Pris: Kontakta oss
Skickas inom 2 veckor

In the field, technicians can test many samples in only a few minutes, and make immediate evaluations on-site. No measuring, set-up, clean-up, or chemical handling are necessary. Lab technicians also use the test strips for pre-test screening-to detect the presence of materials that might interfere with testing. Steps: (mg/L) 0, 0.5, 1.0, 2.0, 4.0, 10.0.

Specifikationer

Enhet:	50 tests
Fotnot :	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Mätområde:	0 - 10 mg/L Cl ₂
Minsta ökningssteg:	Steps: 0, 0.5, 1, 2, 4, 10
Parameter:	Free and Total Chlorine
Plattform:	Teststickor
Test-kit modell:	TS