



Test strips, Nitrate 0 - 50 mg/L and Nitrite 0 - 3 mg/L, 25 tests

Produktnr:

2745425

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Ranges

Nitrate NO₃-N: 0 - 50 mg/L

Nitrite NO₂-N: 0 - 3 mg/L

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Specifikationer

Enhet:	25 tests
Mätområde:	Nitrate 0 - 50 mg/L NO ₃ -N mg/L
Minsta ökningssteg:	Nitrate: 0, 1, 2, 5, 10, 20, 50 mg/L
Parameter:	Nitrate and Nitrite
Plattform:	Teststickor
Smallest increment:	Nitrite: 0, 0.15, 0.3, 1, 1.5, 3 mg/L
Test-kit modell:	Test Strips