



Be Right™



Test strips free chlorine, 0-600 mg/L, 100 tests

Produktnr: 2890200
SEK Pris: Kontakta oss
Tillgänglig

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests

Simplify testing and get results in less than 1 minute, in the field or lab

Eliminate chemical handling and clean-up

Reduce guesswork by screening samples

***Not for use in sea or brackish water**

Specifikationer

Enhet:	100 tests
Fotnot :	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Mätområde:	0 - 600 mg/L Cl ₂
Minsta ökningssteg:	Steps: 0, 25, 50, 100, 200, 400, 600
Parameter:	Free Chlorine, High Range
Plattform:	Teststickor