



Digital titrator assembly

Produktnr:	1690002
SEK Pris:	Kontakta oss
Tillgänglig	

Portable and robust digital Titrator for fast and precise titrations. High resolution (125µL/digit) and high accuracy ($\pm 1\%$). Easy handling. Applications for more than 40 parameters with interchangeable cartridges.

Accurate results - when and where you need them

Fast and Flexible

Versatile

Economical

Specifikationer

Leveransomfattning:	Digital Titrator device only
Tillbehör som ingår:	Nej