



Fast volym, mikropipett, 200 µL

Produktnr:	2946007
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and colour-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Specifikationer

Kapacitet :	200 µL
Noggrannhet:	0,5 %
Volym/förpackningsstorlek:	1 each