



Soxhlet thimbles, cellulose, 33 x 80 mm, 25 pk

Produktnr: 2966700

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Xenosep extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

Xenosep extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

The thimbles are manufactured from high purity alpha cellulose cotton linters which provide excellent mechanical strength, stability, and retention with very low/no background extractables in a disposable, sample holder format.

Specifikationer

Material:	Cellulose
Porstorlek:	10 µm
Volym/förpackningsstorlek:	25 /pk