



## Ice-Pic 20 NTU kalibrerings-/verifieringsmodul

<b>Produktnr:</b>	<b>5225000</b>
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Verify calibration in less than 60 seconds. Available in 1 and 20 NTU modules. Can be used with the 1720C, D and E instruments. No loss of sample flow to complete a verification. No sample cell to recondition after verification.

### **Saves time**

Verification is the fastest available, taking less than one minute.

### **Ease of use**

Compact and lightweight, ready for use at any time.

### **Accuracy**

$\pm 10\%$  or better, and each device is supplied with a Certificate of Accuracy.

### **Repeatability**

Exceptional repeatability performance to be tracked closely.

### **Affordability**

One-time investment, no consumables. One-year warranty.