



Reagens för järn, 0,02-3,00 mg/l Fe

Produktnr:

85428

SEK Pris:

Kontakta oss

Skickas inom 2 veckor

Farligt



Reagens i förseglad förpackning. För 25 ml prov.

FerroVer Iron Reagent: For the determination of total soluble Iron including all ferric (Fe³⁺) and ferrous (Fe²⁺) Iron as well as complexed Iron compounds (Fe-EDTA). Perfect for many applications and a wide range of concentrations. Determines Iron in the presence of Copper.

Specifikationer

Antal tester :	1000
EPA compliant:	Ja
Mätområde:	0,02 - 3,00 mg/L Fe
Metodnamn:	FerroVer
Metodnummer:	8008
Plattform:	Powder Pillows
Volym/förpackningsstorlek:	1000 /pk