



Silvernitrat-patron till Digital Titrator, 1.128 N

Produktnr:

1439701

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Bruksfärdig titreringslösning för bestämning av klorid (silvernitratsmetoden), 10-10000 mg/l Cl⁻

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Mätområde:	100 - 10000 mg/L Cl ⁻
Metodnamn:	Silver Nitrate
Metodnummer:	8207
Parameter:	Chloride
Plattform:	Digital Titrator
Titrant:	1.128 N AgNO ₃