



TitraVer EDTA digital titratorpatron, 0,800 M

Produktnr:

1439901

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Ready-to-use titrant solution for the determination of Calcium and Total Hardness, 100-4000mg/L CaCO₃.

Specifikationer

Antal tester :	40 - 100
Mätområde:	100 - 4000 mg/L CaCO ₃
Metodnamn:	EDTA (CalVer® 2)
Metodnummer:	8204
Parameter:	Hardness, Calcium and Total
Plattform:	Digital Titrator
Titrant:	0.800 M EDTA