



SulfaVer 4 Reagens för sulfat, 2-70 mg/l SO₄

Produktnr:

1206599

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Reagenspulver i tillslutna förpackningar

Specifikationer

Antal tester :	100
Beskrivning:	SulfaVer 4 Sulfate Reagent Powder Pillow, 25 ml
EPA compliant:	Ja
Fotnot :	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Instrument:	DR/4000, DR/3000, DR/2010, DR/2000
Mätområde:	2 - 70 mg/L SO ₄
Metodnamn:	SulfaVer 4
Metodnummer:	8051
Parameter:	Sulfate
Plattform:	Powder Pillows
Reference method:	USEPA 375.4
Volym/förpackningsstorlek:	100 /pk