



ARSENIC TEST KIT

Produktnr: 2822800

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Hach's most economical test kit for determination of arsenic levels.

Simple, economical kit for measuring arsenic levels in the field

Provides fast, reliable results in 20 minutes

Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kapslingstyp:	D
Mätområde:	0 - 500 ppb As
Minsta ökningssteg:	Steps: 0, 35, 75, 175, 1,500, 4,000 ppb
Parameter:	Arsenic
Plattform:	Teststickor
Smallest increment:	Steps: 0, 10, 25, 50, 250, 500 ppb
Test-kit modell:	Test Strip (EZ)‡
Vikt :	1,13 kg