



Silica Standard Solution, 500 µg/L SiO₂, 2.9 L

Produktnr:	2100803
SEK Pris:	Kontakta oss
Tillgänglig	

Silica Standard Solution, 500 µg/L SiO₂ (2.9 L) used in the Series 5000 Silica Analyser.

Specifikationer

Instrument:	5500sc Silica
Koncentration:	500 µg/L
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Parameter:	Silica
Volym/förpackningsstorlek:	2.9 L