



Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Produktnr:	6788302
SEK Pris:	Kontakta oss
Tillgänglig	

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Specifikationer

Produktsats:	5500 sc Silica Analyser, 2/4 Channels
Vad ingår i leveransen:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Vad finns i förpackningen

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions