



Manual flow cell kit for 2100N/AN laboratory turbidimeters

Produktnr: 4744900

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Low-pressure, flow-cell kit allows for a constant optical path.

Obtain colorimetric measurements with Hach spectrophotometers and turbidimeters without individual sample cells. Eliminate errors due to scratched or smudged sample cells. No need for optically matched sample cells. For rapid, convenient processing of large numbers of samples.

Specifikationer

Leveransomfattning:	Flow-cell stand assembly, flow-cell inlet reservoir, reservoir cover, collection drain assembly, flow-cell assembly, interconnecting tubing, flow-cell light shield
Provvolyt:	120 mL purge volume
Vikt :	1,9 kg