



Flow-through cell kit, 230 V

Produktnr: 4745002

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Low-pressure flow-cell kit allows for a constant optical path. Sample flow controlled by a flow valve module.

Obtain colorimetric measurements with Hach spectrophotometers and turbidimeters without individual sample cells. Eliminate errors due to scratched or smudged sample cells. No need for optically matched sample cells. For rapid, convenient processing of large numbers of samples.

Specifikationer

Certifieringar:	TUV/CE
Ingångsspänning:	230 V AC
Leveransomfattning:	Flow-cell stand assembly, flow-cell inlet reservoir, reservoir cover, collection drain assembly, flow-cell assembly, interconnecting tubing, flow-cell light shield
Output Voltage:	12V DC at 1 A
Provvolum:	120 mL purge volume
Vikt :	3,2 kg