



PHC3359 Kombinerad pH-elektrod för mikroprover, fiber, BNC

Produktnr: E16M343
SEK Pris: Kontakta oss
Inget leveransdatum anges

pH electrode with Ag/AgCl reference system. Measuring range 0-12 pH, 0-80°C. Liquid junction: fibre. With BNC plug.

Specifikationer

Garanti :	6 månader
Kompatibla instrument:	MeterLab & sensION+ GLP lab
Kontakt:	BNC Connector
Parameter:	pH
Typ:	Lab, kombinationselektrod