



Differential 1.5" pH sensor, LCP, 5-wire, 7.5m cable

Produktnr:

6028P025

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Differential pH sensor with patented measurement technology using 3 electrodes. Low maintenance, versatile mounting styles, great accuracy and reliability. Chemical resistant body.

Specifikationer

Drift:	0.03 pH per 24 hours, non-cumulative
Elektrodmateriäl:	Glass
Flöde:	max. 3 m/s
Givargänga:	1.5" NPT at both ends
Kabellängd:	7.5 m 25 '
Längd:	180 mm
Material:	LCP
Mätområde:	0 - 14 pH
Monteringsform:	Convertible (immersion or flow-through mount)
Omgivningstemperatur:	-5 - 95 °C
Temperaturgivare:	NTC 300 Ω thermistor for automatic temperature compensation and analyzer temperature readout
Tryckområde:	max. 10.7 bar sensor only (pressure range of mounting hardware to be checked separately)