



pH test kit, 6.5 - 8.5 pH, model 17H

Produktnr: 147008
SEK Pris: Kontakta oss
Inget leveransdatum anges

Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	200
Kapslingstyp:	B
Mätområde:	6.6 - 8.4 pH
Metodnamn:	Color disc/Phenol red
Modell:	17H
Parameter:	pH, mid range
Plattform:	Color Disc
Smallest increment:	0.1 pH
Test-kit modell:	17H