



pH test kit, 7.4 - 9.6 pH, model 17J

Produktnr:

147009

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	200
Kapslingstyp:	B
Mätområde:	7.4 - 9.6 pH
Metodnamn:	Color disc/Thymol blue
Modell:	17J
Parameter:	pH, mid range
Plattform:	Color Disc
Smallest increment:	0.1 pH
Test-kit modell:	17J