



## Test kit, molybdate, model MO-2

**Produktnr:**

**1419301**

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

**Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method**

**Effectively monitors corrosion inhibitors**

---

### Specifikationer

Antal tester :	100
Kapslingstyp:	F
Mätområde:	0 - 8 mg/L Mo
Metodnamn:	Color disc/ Mercaptoacetic acid
Modell:	MO-2
Parameter:	Molybdate - As Molybdenum
Plattform:	Color Disc
Smallest increment:	0.2
	1
Test-kit modell:	MO-2