



## Total Hardness Test Kit, Model 5-EP

**Produktnr:**

**145400**

**Farligt**

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Water hardness is caused almost entirely by calcium and magnesium ion. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With the Model 5-EP test kit, simply count the number of drops required to change the solution color to determine hardness concentration. Utilizes drop count titration with a range of 1-30 gpg. Contains approximately 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Specifikationer

Antal tester :	100
Kapslingstyp:	B
Mätområde:	1 - 30 CaCO <sub>3</sub>
Metodnamn:	Drop count titration/ EDTA-powder pillows
Parameter:	Hardness, Total - As CaCO <sub>3</sub>
Plattform:	Drop Count Titrator
Smallest increment:	1 gpg
Test-kit modell:	5-EP