



Free & Total Chlorine Test Kit (Model CN-70)

Produktnr:

1454200

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	100 Free, 100 Total
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kapslingstyp:	D
Mätområde:	0.02 - 0.68 mg/L Cl ₂
Metodnamn:	Color disc/DPD
Parameter:	Chlorine (Free and Total)
Plattform:	Color Disc
Smallest increment:	0.02
	0.1
Test-kit modell:	CN-70 ²