



## Test kit, total chlorine, model CN-DT, 20-2000 mg/L

**Produktnr:**

**2471100**

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons. Kit utilizes the Hach Digital Titrator and contains 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

**Kit contains all required reagents and apparatus to perform analysis in the field**

---

### Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Kapslingstyp:	F
Mätområde:	20 - 2,000 mg/L Cl <sub>2</sub>
Metodnamn:	Digital Titrator/Thiosulfate
Parameter:	Chlorine, mid and high range - As Total Cl <sub>2</sub>
Plattform:	Digital Titrator
Smallest increment:	0.2 - 5.0
Test-kit modell:	CN-DT