



Test kit, sulfite, model SU-DT

Produktnr:

2063300

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Farligt



Using the Digital Titrator, this kit contains 100 tests with a range of 10-800 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	F
Mätområde:	4 - 800 mg/L SO_3^{2-}
Metodnamn:	Digital Titrator/Iodometric
Minsta ökningssteg:	0.4-4.0**
Parameter:	Sulfite - As SO_3^{2-}
Plattform:	Digital Titrator
Test-kit modell:	SU-DT