



Test kit, alkalinity, model AL-DT

Produktnr:

2063700

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Kapslingstyp:	F
Mätområde:	10 - 4,000 mg/L CaCO ₃
Metodnamn:	Digital Titrator/Sulfuric acid
Minsta ökningssteg:	0.1-10
Parameter:	Alkalinity - (P) & Total (MO)
	as CaCO ₃
Plattform:	Digital Titrator
Test-kit modell:	AL-DT