



Test kit, alkalinity, model AL-TA

Produktnr:

2314500

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kapslingstyp:	B
Mätområde:	23 - 495 gpg & 385 - 8500 mg/L
Metodnamn:	Drop count titration/Sulfuric acid
Parameter:	Alkalinity - (P) & Total (MO) As CaCO ₃
Plattform:	Drop Count Titration
Smallest increment:	23 gpg
	385
Test-kit modell:	AL-TA