



Test kit, ozone, model OZ-2

Produktnr:

2064400

SEK Pris:

Kontakta oss

Skickas inom 2 veckor

Farligt



Hach's most economical test kit for ozone measurement utilizes the simple, low-cost color disc method.

Hach's most economical test kit for ozone measurement

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	100
Kapslingstyp:	B
Mätområde:	0 - 2.2 mg/L O ₃
Metodnamn:	Color disc/DPD
Modell:	OZ-2
Parameter:	Ozone
Plattform:	Color Disc
Smallest increment:	0.05
Standardmätområde:	0 - 2.2 mg/L O ₃
Test-kit modell:	OZ-2