



Nitrogen, Nitrite Colour Disc Test Kit, Model NI-6

Produktnr:

224000

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Utilising the simple colour disc method/ferrous sulfate, this kit measures Nitrogen, Nitrite as NO_2^- and contains 100 tests. Range (mg/L): 0-8/0-1600.

Simple, low-cost colour disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester : 100

Kapslingstyp: D

Mätområde: Single Range

Metodnamn: Colour disc/

Parameter: Ferrous sulphate

Plattform: Nitrogen, Nitrite

Smallest increment: Color Disc

2

40

Test-kit modell: NI-6

Vad ingår i leveransen: Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.

Vad finns i förpackningen

Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.