



## Cylinder, graduated, double metric, 1 L

|                           |                |
|---------------------------|----------------|
| <b>Produktnr:</b>         | <b>2928305</b> |
| SEK Pris:                 | Kontakta oss   |
| Inget leveransdatum anges |                |

These glass pyrex@ graduated cylinders feature two sets of graduations for reading either up or down Calibrated "to deliver"

---

### Specifikationer

|                            |               |
|----------------------------|---------------|
| Certifieringar :           | Class A       |
| Höjd :                     | 465 mm        |
| Kapacitet :                | 1000 mL       |
| Material:                  | Glass         |
| Mätmetod:                  | To Deliver    |
| Typ:                       | Double Metric |
| Volym/förpackningsstorlek: | 1 each        |