



## Beaker, polypropylene, low form, 50 mL, 12/pk

<b>Produktnr:</b>	<b>108071</b>
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Specifikationer

Diameter, yttre:	51.1 mm with spout
Enhet:	12 /pk
Kapacitet :	50 mL
Material:	Polypropylene, Low Form