



CEL/890 Advanced portable laboratory

Produktnr:

2994000

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



Offering on-site results in an easy-to-use, rugged, portable lab, Hach's Oil & Gas Water Quality Laboratory covers parameters critical to effective source water and base water analysis. Depending on water composition, this kit may also be used to support analysis of produced water, flowback water, water treatment, and enhanced oil recovery applications.

Specifikationer

Applikation:

Fracing

Vad ingår i leveransen:

This laboratory includes: DR 900 Colorimeter, HQ40d Multi-Parameter Digital Meter, Digital Titrator, all required reagents for 12 critical parameters, apparatus and accessories, and two convenient, rugged carrying cases.

Vad finns i förpackningen

This laboratory includes: DR 900 Colorimeter, HQ40d Multi-Parameter Digital Meter, Digital Titrator, all required reagents for 12 critical parameters, apparatus and accessories, and two convenient, rugged carrying cases.