



ORP platinum electrode 5362, gel electrolyte, 3 diaphragms, 0-80 °C, 6 bar, S8 screw head

Produktnr: LZU5362.97.0002
SEK Pris: Kontakta oss
Inget leveransdatum anges

ORP platinum electrode 5362 with gel electrolyte and S8 screw head. ± 2000 mV, 0 - 80 °C, up to 6 bar. 3 ceramic diaphragms.

Suitable for swimming pools,

drinking water, gas scrubbing towers, boiler water, detoxification, wwtp outlets etc.

Specifikationer

Elektrolyt:	Gel
Garanti :	6 månader
Kabelanslutning :	S8
Längd:	120 mm
Material:	Sensor: glass/Pt
Mätområde Redox:	± 2000 mV
Membran:	3 ceramic
Reference system:	Encapsulated Ag/AgCl crystals
Temperaturområde :	0 - 80 °C
Tryckområde:	upp till 6 bar