



## Conductivity cell 5396, k=1, 80 °C, S8 connector, PG 13.5

|                           |                        |
|---------------------------|------------------------|
| <b>Produktnr:</b>         | <b>LZU5396.97.0002</b> |
| SEK Pris:                 | Kontakta oss           |
| Inget leveransdatum anges |                        |

Conductivity cell 5396, in glass and platinum, k=1, 1...200000 µS/cm, 0...80 °C. S8 connector, PG 13.5 thread. (Without A.T.C.)

---

### Specifikationer

|                         |                          |
|-------------------------|--------------------------|
| Cellkonstant k:         | 1 cm <sup>-1</sup>       |
| Kabelanslutning :       | S8                       |
| Längd:                  | 110 mm                   |
| Material:               | Measuring Cell: Glass/Pt |
| Mätområde:              | 1 - 200.000 mg/L (µS/cm) |
| Temperaturkompensation: | 0 - 80 °C                |