



Conductivity cell 5398 with PT1000, k=0.3, 80 °C, MP-5, PG 13.5

Produktnr: LZU5398.97.0002

SEK Pris: Kontakta oss

Inget leveransdatum anges

Conductivity cell 5398 + PT1000 A.T.C., titanium, k=0.3, 0.5...80000 µS/cm, 0...80 °C. MP-5 connector, PG 13.5 thread.

Specifikationer

Diameter:	12 mm
Garanti :	6 månader
Kabelanslutning :	MP-5
Längd:	110 mm
Material:	Measuring Cell: Titanio
Mätområde:	0.50 - 80.000 mg/L (µS/cm)
Temperaturkompensation:	0 - 80 °C
Temperaturområde :	Pt 1000
Vad ingår i leveransen:	Includes: sensor with cable and manual

Vad finns i förpackningen

Includes: sensor with cable and manual