



EVITA OXY 1100 Membrane 25µm for Dissolved Oxygen sensor, 1 pc

Produktnr:

085G0021

SEK Pris:
Tillgänglig

Kontakta oss

Recommended for measuring range 0.002-2mg/L DO.

Specifikationer

Avfall:	OXY 1100 sensor must be disposed as ordinary industrial waste.
Mätområde:	0,002 - 2 mg/L
Responstid :	7 s
Tjocklek, membran:	25 µm
Volym/förpackningsstorlek:	pk/1