



Be Right™



APA STD 2, LR HDNS 2000 UG/L 1L

Produktnr:	2696353
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

APA6000 Low Range Hardness Standard 2, 5 mg/L, 1L

Specifikationer

Instrument:	APA6000
Koncentration:	5 mg/L
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Parameter:	Hardness
Volym/förpackningsstorlek:	1 L