



## APA6000 Alkalinity Cleaning Solution

<b>Produktnr:</b>	<b>2697053</b>
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

The APA6000 Alkalinity Process analyzer uses five reagents each of which is supplied in a 1L bottle designed to provide reagents for 30 days of continuous use.

**APA Alkalinity Reagent 1 (2826153) is an acid titrant that reacts with the alkalinity in the sample.**

**APA Alkalinity Reagent 2 (2696653) is an indicator used to determine the endpoints of the titration.**

**APA Alkalinity Standards 1 and 2 (6001100) are used together to perform a multiple point calibration.**

**The fifth reagent, APA Alkalinity Cleaning Solution (269703), cleans the system during the Prime and Instrument Clean functions.**

---

### Specifikationer

Instrument:	APA6000
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Mätområde:	N/A
Parameter:	Alkalinity
Volym/förpackningsstorlek:	1 L