



## PROBE, DISSOLVED OXYGEN, 3M CABL

Produktnr: 5197003

### UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Clark-Type Polarographic Dissolved Oxygen probe with 3 meter cable and 5-pin circular connector compatible with Hach Sension meters.

#### Accurate

You don't just need an answer, but the right one, and fast.

#### Easy to Use

Now, more than ever, you're looking for new ways to stretch your dollar and resources. Skip studying complicated user manuals and spend your time measuring, thanks to Hach's innovative menu designs and simple keypads for truly intuitive operation.

#### Flexible

Need to measure in both the lab and field? Hach offers efficiency with reliable, high-performance systems that easily move with you from the benchtop directly to the testing site. Measure the parameters you need today with the assurance you can expand them tomorrow.

#### Rugged and Durable

Rely on Hach's durable measuring systems designed to withstand the rigors of outdoor testing and frequent use. Choose from the industry's toughest probes or non-glass electrodes for rugged, waterproof systems that deliver day after day.

---

### Specifikationer

Elektrolyt:	2759126
Garanti :	6 månader
Givartyp:	Platinum cathode, silver anode with replaceable PTFE membrane module (item no. 5197300)
Kabellängd:	3 m
Längd:	149 mm
Material:	Sensor Body: Epoxy
Mätområde:	0 - 20 mg/L
Temperaturområde :	0 - 50 °C

---

### Obligatoriska tillbehör

- Påfyllningslösning, intern, polarografiskt DO, 50 mL (Item 2759126)