



## CDC241-9 Konduktivitetscell

<b>Produktnr:</b>	<b>E61M008</b>
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Conductivity cell for continuous measurements, titration, viscous media.

### **The right conductivity cell for your application.**

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.

---

### **Specifikationer**

Conductivity Cell:	CDC241-9
Elektrod:	2-pole, non-platinised
Garanti :	12 månader
Specialfunktioner:	Conductivity cell for continuous measurements, titration, viscous media.