



GK2401C Kombinerad pH-elektrod, röd stav, poröst stiftdiaphragma

Produktnr: 945-252
SEK Pris: Kontakta oss
Tillgänglig

Long-term accuracy and reproducibility with Red Rod technology

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100°C. Liquid junction: porous pin.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

Specifikationer

Dimension (D x L):	9.5 mm x 103 mm
Elektrolyt:	Refillable: Saturated KCl
Garanti :	12 månader
Givar koppling:	Porous pin
Givartyp:	Glass
Kabellängd:	1 m
Längd:	103 mm
Material:	Glass
Mätområde:	0 - 12 pH
Temperaturområde :	-10 - 100 °C
Termistor:	none

Obligatoriska tillbehör

- Påfyllningslösning, referens, mättad KCl (KCl.L), 100 mL (Item S21M002)